



## Product Catalog

2007



For All Of Your Beauty Desires

- *Electrolysis*
- *Skin Care*
- *Vascular*
- *Accessories*

## Table of Contents:

### Electrology Systems **4**

- Elegance
- Essence
- Nova 2000

### Skin Care Systems **6**

- Jeunesse
- Hydrophoresis
- Swiss Spa System

### Vascular System **9**

- Skin System

### Skin Care Lotion and Gels **10**

### Accessories **11**

- Electrology **11**
- Skin Care **13**
- Microcurrent **19**
- Vascular **20**
- Miscellaneous **21**

## CD/DVD's/Brochures



Jeunesse Glove Procedure DVD  
Order# DVD-2006



Jeunesse Waiting Room DVD  
Order# DVD-2006WR



Telangiectasia DVD  
For use with the Skin System  
Order# DVD-2007



Microcurrent History, Studies,  
Testimonials, Before/After  
Pictures CD  
For use with the Skin System  
Order# CD-2007

# Miscellaneous items



The Blend Method by Michael Bono  
Learn the concept of blend from theory to application.  
Order# BO-0001

Telangiectasia by Michael Bono  
Learn the removal of telangiectasia (broken capillaries)

Order# BO-0002




Telangiectasia Needles  
For use with the Skin System  
Order# NE-00K3



Jeunesse  
Waiting  
Room  
Brochure

Order#  
BRO-2006



Patient History Card  
Record your client's medical history and treatments. Back of card provides area for listing treatment settings and face and body figure for noting treatment areas and problems. Additional back pages can be purchased as needed for client's file.

Order# 8900-140



**3555 Airway Drive  
Suite 307  
Reno, Nevada 89511  
Toll Free 1-800-334-7126  
Fax 775-332-3852  
[www.clareblend.com](http://www.clareblend.com)  
[staff@clareblend.com](mailto:staff@clareblend.com)**

## Electrology



The **Elegance** is Clareblend's newest epilator. Digital technology provides precision current application. Pre-programmed settings for blend are available for every hair type. With one touch of a button the displays show units of lye, DC, RF, and timer settings. Manual adjustments are provided for each modality. Preference buttons allow you to store your most used settings. Choose your Preference button and settings are recalled on displays and system is ready for operation. Most importantly the technician can manually adjust any of the pre-programmed settings. The machine is under your control.



The **Essence**, low intensity epilator, offers single needle manual operation of blend, galvanic, and self-seeking Autoblend. Phoresis modality included. New Audio treatment timer. The Essence is small and lightweight, making it perfect for the Electrologist who needs mobility.

## Miscellaneous Accessories



Extra Large Case  
Order# 8900-132

Large Case  
Order# 8900-130

Small Case  
Order# 8900-131



Equipment Cart with Power  
strip Built In

Three Melamine shelves  
One drawer for accessories  
Top shelf 16"x 22", with  
rounded corners  
Power bar with 6 electrical  
outlets  
Safety circuit breaker  
Lamp support on each side of  
table  
Partial assembly required  
Dimensions: 22" L x 18.5" W x  
38" H  
Weight 36lbs

Order# 9503-001

## Accessories continued



Microcurrent Extender Cords

Order# 2006-001

## Acupuncture Accessories



Microcurrent Acupuncture cord

Order# 2004-002



RF Acupuncture cord

Order# 2004-004



*Multiple-needle epilator with medical grade armature.*

Clareblend is the only manufacturer to offer the BLEND modality on a multiple-needle epilator, Single DC and Single RF current included.

The **Nova 2000®** is fast and effective. Clareblend's patented electronics has significantly reduced multi-needle epilation time. The client will experience a reduced treatment time, more comfort and quicker satisfaction. For the technician's comfort, the indicator lights, timer, DC/Blend controls and DC meter have been mounted on the end of the armature. The medical grade armature reduces the problem of probes pulling out of the follicles.

Each of the 6 channels, or probes, are individually controlled. An automatic two second delay gives the technician the time needed for insertion. The slow rise gives the customer a comfortable treatment. Simplistic innovative design allows the technician to master the Nova 2000® very quickly.

## Skin Care



The **Jeunesse** is a microcurrent system for skin rejuvenation and anti-aging.

**Benefits:**

- Increase collagen and elastin
- Increase ATP production
- Increase skins firmness
- Diminish fine lines and wrinkles

**Modalities include:**

- Fixed Frequency - HI 0-1000 Hz
- Mixed Frequency - 1-50 Hz
- Micro current 0 to 1000 micro amps
- Random modulation
- 3 Pre-sets: Face, Mouth/ Eyes, Body
- Session Timer
- Conductivity Tester
- Option:
  - LED Red or Blue Light Therapy



The **Swiss Spa System**

is a gentle low-level galvanic system for home or professional use. The current combined with the Spa Gels will reduce fine lines and wrinkles, increase production of collagen and elastin providing the look of firm and more youthful skin.

## Accessories continued



Footswitch Extensions  
Extend the reach of your footswitches by 2 feet.

Order# 9501-050

## Microcurrent Accessories



Microcurrent Handles  
For facial rejuvenation procedure

Order# 2004-003



Microcurrent treatment cord

Order# 2004-001

# Accessories continued



Silent Foot Switch

Order# 8900-014



Round Foot Switch

Order# 8900-015



Footswitch "Y" Adapter

Allows you to use 1 footswitch for blend modality and have individual control of each meter.

Order# 9501-051



The Hydrophoresis is a low-level galvanic treatment system for all skin types. It boosts the penetration of product into the skin, maximizing its performance. The gentle galvanic current and massaging action of the spheres increases circulation, helps to hydrate and soothe the skin leaving it with a natural healthy glow and radiance. It also prepares the skin by softening the keratin and sebum, which assists in the removal of blackheads and congestion.

# Vascular System



The **Skin System** unit treats vascular and skin blemishes.

Modalities include:

- Pre-programmed settings for small, medium & large Telangiectasia
- Hi & Low output RF current
- DC current
- Cata (+) and Ana (-) Phoresis
- 1 foot pedal operation

# Accessories continued



Europe Power Cord

Order# 8900-071



UK Power Cord

Order# 8900-072



Australia Power Cord


Order # 8900-073

## Accessories continued



Probe cord tips  
End Bibs - cone shaped

Order# 8900-170



Probe cord tips  
End Bibs - short tips

Order# 8900-171



U.S. Power Cord

Order# 8900-070

## Gels & Lotions



- **Clareblend DMAE/Hyaluronic Gel**
- Conductive gel for use with Jeunesse and Clareblend galvanic systems.



### Dr. Yao

Herbal Acne Face wash and Pimple cream. When used with the galvanic and microcurrent it may help to control the sebum and acne breakouts.

Dr. Yao's cosmetics unite the theory of Traditional Chinese Medicine with modern achievements of science making skin soft, luminous and beautiful.

## Gels & Lotions



The **"Pre-treat"** lotion increases collagen and elastin; firms skin, fights free radicals, improves skin tone/texture, softens, soothes and enhances bio balance. Increases ATP production.

The **"Treatment"** lotion revitalizes skin, aides in reducing fine lines and wrinkles, fights free radicals, improves lipid barrier, increases moisture and radiance, balances pH.

2 each Pre-Treat Gel (negative polarity)  
2 each Treatment Gel (positive polarity)

Used with the Swiss Spa System, Hydrophoresis, or Jeunesse.



### SesDerma Ioson DMAE Gel

Reaffirming gel for use with the Jeunesse, Swiss Spa System, and Galvanic Systems.. Creates a lifting effect and fights flaccidity, increasing muscular tone. It attenuates and softens wrinkles, reaffirming and improving the elasticity of the skin.

## Accessories continued



Phoresis Wand Electrode

Order# 9001-240



Accessory Cords

Order# 8900-050



Rubber Holders

Order# 8900-060

## Accessories continued



Phoresis Facial Roller

Order# 8900-030



Phoresis Body Roller

Order# 8900-040



Phoresis Sphere Electrode

Order# 9001-250

## Electrology Accessories



Probe Cord for Needles

1 pc K-Shank 8900-001  
2 pc F-Shank 8900-002  
Taper 8900-007



Multineedle Probe Cord  
Nova 2000™ 24" cord

Order# 9506-091



Ground Electrode-Small  
Aluminum


Order# 8900-020

## Accessories continued



Ground Electrode-  
Large  
4" x 1" diameter  
Medical grade  
Stainless Steel

Order# 8900-025



Flat Ground  
Electrode-Small  
2" x 4"  
Stainless Steel

Order# 8900-022



Flat Ground  
Electrode-Large  
4" x 4"  
Stainless Steel

Order# 9601-0GP

## Skin Care Accessories



Conductive Gloves

Order# 2006-05G



Snap Cords for  
Conductive Gloves

Order# 2006-002



LED Light Therapy

Order# LED-2006B  
Blue Light

Order# LED-2006R  
Red Light

Option to add to  
Jeunesse or  
Elegance Systems